

Line-of-sight port ●  
<sup>3</sup>He dump ●

Universal housing ●

# HelioxVT™

## Why choose HelioxVT?

The HelioxVT range of single shot <sup>3</sup>He systems allows users to access temperatures below 300 mK for extended periods.

A fully configured MercuryiTC provides total control of the HelioxVT, automating cool down from room to base temperature and simplifying integration into your measurement setup via a range of standard communication interfaces.

For more specific experimental requirements, we can offer tailored <sup>3</sup>He systems designed to meet your needs.

## Precise control of magnetic field and temperature


The HelioxVT is designed to operate safely in integrated into cryo-magnet systems – allowing access to the lowest temperatures and the highest fields.

Sorption pump ●

1 K condensing stage ●

IVC flange ●

<sup>3</sup>He pot ●



## Features

- Achieves < 300 mK for > 40 hours and 50  $\mu$ W of cooling power at 350 mK for > 6 hours
- Fast turn around time for sample exchange
- **HelioxVT** uses a cold gas environment with a 50 mm access, therefore no liquid helium
- The sample temperature range of a new or existing VTI can be extended below 300 mK
- No liquid helium in the sample plane making the **HelioxVT** ideal for neutron or X-ray scattering experiments
- Compatible with 50 mm diameter VTIs
- 1 K pot free design – no additional room temperature pumps make this a simple, self-contained solution.

## Magnetic Field Configuration

Magnetic field requirement	Configuration	Benefits
Up to 14 T	HelioxVT with TeslatronPT Cryofree superconducting magnet system	<ul style="list-style-type: none"><li>- No requirements for liquid cryogen (or accompanying infrastructure)</li><li>- Complete turn-key solution for material characterisation</li></ul>
Up to 21 T	HelioxVT Integra liquid helium cryostat	<ul style="list-style-type: none"><li>- Higher magnetic fields than a Cryofree system</li><li>- IntegraAC cryostat offers almost zero liquid helium consumption (at 4.2 K)</li></ul>

## Key Specifications

Base temperature	$\leq$ 300 mK for 40 hrs with no applied heat load
Cooling power	< 350 mK for 6 hrs with 50 $\mu$ W applied
Temperature range	300 mK to 300 K
Temperature stability	$\pm$ 3 mK below 1.2K ; $\pm$ 0.1 K above 1.2 K
HelioxVT Sample space	43 mm diameter

Visit [nanoscience.oxinst.com](https://nanoscience.oxinst.com) or email [nanoscience@oxinst.com](mailto:nanoscience@oxinst.com)